

**HEINE** QUALITY  
MADE IN GERMANY

## ∴ **HEINE BETA® 400 F.O. Otoscope**

The highest magnification for the best diagnosis.

NEW!

**4.2x**

MAGNIFICATION



DEPTH OF FIELD



ONE-HAND  
OPERATION



# THE HEINE BETA® 400 FIBER OPTIC OTOSCOPE. 4.2x MAGNIFICATION. EXCEPTIONAL DEPTH OF FIELD. ONE HAND OPERATION.

The HEINE BETA400 F.O. Otoscope is the most advanced high quality otoscope making more details visible, while maintaining familiar handling and functionality. With a previously unmatched magnification factor of 4.2x and multi-coated precision optics for sharp spatial images, it truly provides for the best diagnosis. Combined with a large field of view the ear canal and tympanum are clearly visible at a glance. The unique opening mechanism (patent pending) can be operated with one hand and flips up the optical system. Application of instruments through the device becomes easily possible.



## 4.2x

### 4.2x MAGNIFICATION

Unrivalled: More detail of anatomical structures and easy recognition of foreign bodies.



### ONE HAND OPERATION

Easy thumb control: Optical system to be flipped open and out of the way, enabling instrumentation and ease of cleaning.



### EXCEPTIONAL DEPTH OF FIELD

Precise image: Provides sharp spatial images in a convenient working range.

BETA400 F.O. Otoscope **without handle** [01]

**2.5V: B-001.11.400**

**3.5V: B-002.11.400**

BETA400 F.O. Otoscope with 1 set (4 pcs.) of reusable tips and 5 of both AllSpec disposable tips 2.5 and 4 mm Ø, spare bulb, **BETA battery handle 2.5V**, hard case [02]

**B-143.10.118**

**Set:** BETA400 F.O. Otoscope with 1 set (4 pcs.) of reusable tips and 5 of both AllSpec disposable tips 2.5 and 4 mm Ø, **BETA200 Ophthalmoscope head**, 1 spare bulb each, 1x BETAL rechargeable handle 3.5V, NT300 charger, hard case

**A-153.29.420**

**Set:** BETA400 F.O. Otoscope with 5 of both AllSpec disposable tips 2.5 and 4 mm Ø, **DELTA20 Plus Dermatoscope head**, spare bulb for BETA400, 1x BETAL rechargeable handle 3.5V, NT300 charger, zipper case

**A-160.29.420**

Further sets on request.

