

Freedom to Care More

You want to give your patients the very best care and spare them from unnecessary discomfort. 3M™ Kind Removal Silicone Tape is a smarter medical tape with the ideal blend of properties to help you care for your patients who have fragile or at-risk skin.

Care without
Compromise



In recent testing,
9 out of 10 nurses preferred
3M™ Kind Removal Silicone Tape
over their current gentle tape.*

Security Offers reliable yet pliable fixation, remains in place until you decide otherwise.

Comfort Removes cleanly, without disrupting fragile skin layers or causing patients any undue pain.


Ease of Use Can be repositioned and neatly torn by hand.



One of a Kind

Unique characteristics make 3M™ Kind Removal Silicone tape unlike any medical tape you've used before. The silicone adhesive is softer and more conformable than traditional adhesives, and works differently with skin and hair. The tape can be removed with

minimal stripping of skin cells or pulling of hair. And unlike traditional tape that increases adhesion over time, 3M silicone tape delivers constant adhesive strength for as long as it's left in place.




Skin Cell Removal Visualization

Proven Comfort

Research confirms that upon removal, 3M™ Kind Removal Silicone Tape causes minimal epidermal cell stripping and less pulling of hair.

Patients will feel the difference the moment it goes on... and comes off.



3M™ Kind Removal Silicone Tape

Traditional (Acrylate) Soft Cloth Tape

Traditional (Acrylate) Paper Tape

Images show skin cells and hair (in red) that were stripped on removal of three tapes. Representative sample from one individual.

Ordering Information

Catalog No.	NHRC/NDC No.	Size	Rolls/Box	Boxes/Case	HCPCS Code**
Single Use Rolls					
2770S-1	8333-2770-11	1 in. x 1.5 yd. (2,5cm x 1,3m)	100	5	A4452
2770S-2	8333-2770-12	2 in. x 1.5 yd. (5,0cm x 1,3m)	50	5	A4452
Standard Rolls					
2770-1	8333-2770-01	1 in. x 5.5 yd. (2,5cm x 5,0m)	12	10	A4452
2770-2	8333-2770-02	2 in. x 5.5 yd. (5,0cm x 5,0m)	6	10	A4452

* From pre-market clinical evaluations.

** Disclaimer: HCPCS codes have been provided to assist you in the preparation of insurance claims. Please note, however, that the reimbursement information provided by 3M Health Care and its representatives is intended to provide general information relevant to coverage and coding for 3M products. Insurers' reimbursement policies can vary and the use of the codes discussed here does not guarantee that an insurer will cover or pay at any particular level. Health care providers should exercise independent clinical judgment in choosing the codes which most accurately describe the products provided.



Skin & Wound Care Division 3M Health Care

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
USA
1-800-228-3957
www.3M.com/healthcare

3M Canada

P.O. Box 5757
London, Ontario N6A 4T1
Canada
1-800-364-3577
www.3M.com/ca/healthcare



Try a sample today!

Request a sample card at www.3M.com/TryBlue, or contact your local 3M sales representative or the 3M Healthcare Customer Helpline at **1-800-228-3957** for more information.