

HEINE K180[®] Ophthalmoscope Otoscope set

Features

- Corneal and iris reflexes are eliminated.
- XHL Xenon Halogen Technology.
- Bright, white light. Clear retinal image.
- Range of lenses: 27 lenses from -35 D to +40 D.
- Brilliant image even at high corrections.
- Ergonomical shape.
- The instrument fits the orbita comfortably.
- Soft orbital rest.
- Steadies the instrument and protects your glasses.
- Choice of 5 apertures.
- For both large and small pupils.
- Choice of fixation star or cobalt blue filter for corneal examination.
- Dust resistant housing.



SAVE UP TO 25%

[Click here to see on line](#)

medent
MEDICAL

HEAD OFFICE:
PH 04 527 7702 FAX 04 527 7499
PH OUTSIDE WELLINGTON: 0800 26 333 68
TOLL FREE FAX: 0800 52 774 99

HEINE QUALITY
MADE IN GERMANY



HEINE K180[®] Ophthalmoscope Otoscope set



Code	Description	Was	Now	Save
A-279.10.118B	BETA BATTERY HANDLE	\$1099	\$887	\$212



Code	Description	Was	Now	Save
A-279.27.388B	BETA4 USB RECHARGEABLE HANDLE WITH USB CORD AND PLUG-IN POWER SUPPLY	\$1625	\$1,219	\$406



Code	Description	Was	Now	Save
A-279.23.420B	BETA4 NT RECHARGEABLE HANDLE WITH NT4 TABLE CHARGER	\$1832	\$1,490	\$342