



# Orthopaedic & Disability Products



## Product Catalogue

2019-20 | Edition 16



[Online Sales & Ordering](#)



[Email Enquiries](#)

## About Medent Medical

After 30 years supplying the medical profession, Medent Medical remains 100% kiwi owned and 100% committed to delivering to you the best of medical, surgical, pharmaceutical and equipment supplies whenever you need them.

## BISMOS - The innovative inventory management tool

We wanted to make it easier for you to order, store and manage your medical supplies. We came up with BISMOS. A complete stock management and ordering solution that takes the pain out of supply management.

With the help of BISMOS you can:

- ✓ Automatically generate your next order by scanning products in your storeroom
- ✓ Easily check stock levels
- ✓ Improve practice productivity
- ✓ Cut medical supply ordering in half
- ✓ Reduce out of stock items
- ✓ And much more

[Discover how BISMOS can work for you](#)

## Providing you the best of medical, surgical & pharmaceutical supplies

With a constantly expanding inventory with thousands of products to choose from, including all of the major brands, when it comes to [medical supplies](#) think of Medent Medical as your 'go to' company. We have been supplying the medical profession with the best of medical, surgical, pharmaceutical and equipment supplies for more than 30 years. That is 30 years of working for you to simplify your ordering needs. You simply won't find a faster, more efficient way to manage your inventory.

Ask how we can help you:


- ✓ Optimise your ordering process
- ✓ Streamline inventory replacement
- ✓ Facilitate easy online ordering
- ✓ Enhance your productivity
- ✓ Improve ordering efficiency
- ✓ And save yourself both time and money

For more choice, better service and expert advice on your medical supplies, [contact us today](#)

# BISMOS AUTOMATED ORDERING AND STORAGE SOLUTIONS



- ✓ Save your time spent ordering supplies
- ✓ Reduce order error rate
- ✓ Easy checking of what stock you have
- ✓ Flexible modular system



## TABLE OF CONTENTS

<b>CATEGORY: BANDAGES-CAST</b> .....	1
BSN Delta Cast Elite .....	1
BSN Delta Cast Soft .....	2
BSN Delta Dry Cast .....	2
BSN Dynacast Prelude .....	3
BSN Gypsona* BP Plaster of Paris Bandages .....	3
3M Scotchcast Plus Casting Tape .....	4
3M Scotchcast Softcast Tape .....	4
Primenet Polyester Stockinette .....	5
<b>CATEGORY: ORTHOPAEDIC PADDING</b> .....	6
BSN Soffban Orthopaedic Padding .....	6
Primepad Polyester Undercast Padding (Soffban Alt) .....	6
<b>CATEGORY: CASTING PROTECTOR</b> .....	7
Aquashield Waterproof Cast Protector .....	7
<b>CATEGORY: BRACE/SPLINT/SUPPORT</b> .....	8
Walker Moon Boot .....	8
Ankle Stabiliser .....	8
AFO-Deluxe-Allcare .....	8
BSN Cast Shoe Adult .....	9
Arm Sling - Immobiliser with Waist Strap .....	9
Elastic Knee Support/Stabiliser .....	9
Sling Arm with Strap Support .....	9
Mesh Arm Sling .....	10
Neoprene Wrist Stabiliser .....	10
Universal Wrist Splint .....	10
Thumb brace .....	10
Finger Splint Zimmer (Aluminium) .....	11
Collar Neck Hard Adjustable Adult .....	11
Finger Splints Stax .....	11
Cervical Collars .....	11
Cardboard Splint .....	12
<b>CATEGORY: MOBILITY</b> .....	13
Crutches Adult Double Adjustable Forearm (Pair) .....	13
Stoppers for Crutches (Pair) .....	13

## CATEGORY: BANDAGES-CAST

### BSN Delta Cast Elite

#### Description

BSN Delta Cast Elite is a polyester cast tape designed with ease of application and patient comfort in mind. The patented fibreglass-free casting material combines strength and rigidity with soft edges for enhanced patient comfort.

Provides rigid casting options for primary and secondary immobilisation, where strength and radiolucency are required.

Patented polyester substrate provides exceptional moulding, and the unique resin system ensures maximum lamination and end laydown for long lasting casts.

#### Features and Benefits

- ✓ Smooth cast with soft edges
- ✓ Fibreglass free
- ✓ Clean, easy to apply and remove
- ✓ Highly conformable
- ✓ Strong yet lightweight
- ✓ Durable to reduce cast breakdown
- ✓ Smooth finish
- ✓ Consistent set times
- ✓ Enhanced radiograph viewing



Medent Code	Mft Code	Size	Unit of Sale	Medent Code	Mft Code	Size	Unit of Sale
10063	72279-01	White 7.5cm x 3.6m	Box 10	10062	72279-03	White 5cm x 3.6m	Box 10
10004	72278-02	Black 10cm x 3.6m	Box 10	10064	72279-02	White 10cm x 3.6m	Box 10
10012	72276-00	Blue 5.0cm x 3.6m	Box 10	10065	72279-03	White 12.5cm x 3.6m	Box 10
10013	72276-01	Blue 7.5cm x 3.6m	Box 10	10032	73455-00	Pink 5cm x 3.6m	Box 10
10053	72275-01	Red 7.5cm x 3.6m	Box 10	10033	73455-01	Pink 7.5cm x 3.6m	Box 10
10054	72275-02	Red 10.0cm x 3.6m	Box 10	10034	73455-02	Pink 10cm x 3.6m	Box 10
10043	72277-01	Purple 7.5cm x 3.6m	Box 10	10022	73454-00	Green 5cm x 3.6m	Box 10
10052	72275-00	Red 5cm x 3.6m	Box 10	10023	73454-01	Green 7.5cm x 3.6m	Box 10
10002	72278-00	Black 5cm x 3.6m	Box 10	10024	73454-02	Green 10cm x 3.6m	Box 10
10003	72278-01	Black 7.5cm x 3.6m	Box 10	10072	72273-22	Camo 5cm x 3.6m	Box 10
10014	72276-02	Blue 10cm x 3.6m	Box 10	10073	72273-06	Camo 7.5cm x 3.6m	Box 10
10042	72277-00	Purple 5cm x 3.6m	Box 10	10082	72273-17	Pastel 5cm x 3.6m	Box 10
10044	72277-02	Purple 10cm x 3.6m	Box 10	10083	72273-07	Pastel 7.5cm x 3.6m	Box 10

## BSN Delta Cast Soft

### Description

BSN Delta-Cast Soft is fibreglass-free and was designed to provide variable support to fractures and soft tissue injuries. Its semi-rigidity allows for more flexibility and mobility, thus helping to reduce immobilisation related problems such as atrophy or joint stiffness. As a result, the overall healing process tends to be shortened.



### Features and Benefits

- ✓ Earlier patient mobility as cast will allow a limited range of movement
- ✓ Improved circulation, faster healing, shorter rehabilitation period
- ✓ Help avoid or minimise muscle atrophy and other immobilisation related issues
- ✓ Lightweight, breathable material offers high wearing comfort
- ✓ Smooth surface and edges for a snag-free finish
- ✓ Available in a wide range of attractive colours

Medent Code	Mft Code	Size	Unit of Sale
10102	72708-08	Red 5.0cm x 3.6m	Box 10
10092	72708-05	Blue 5.0cm x 3.6m	Box 10
10103	72708-18	Red 7.5.0cm x 3.6m	Box 10
10094	72708-25	Blue 10.0cm x 3.6m	Box 10
10111	72708-02	White 2.5cm x 3.6m	Box 10
10112	72708-12	White 5cm x 3.6m	Box 10

## BSN Delta Dry Cast

### Description

The unique patented technology in BSN Delta-Dry® has created a water resistant cast padding which allows easy application by the clinician and affords the patient all the lifestyle benefits of a cast they can get wet.

Delta-Dry® is a water resistant padding while Delta-Dry® Softliner is a water resistant padded stockinette and can be used without additional padding.

### Features and Benefits

- ✓ Better long term water resistance
- ✓ Offers more under cast protection
- ✓ Eliminates the need for stockinette
- ✓ Allows patients to shower, swim and bathe like normal
- ✓ No special tools are needed for application or removal
- ✓ Conformable gives patient increased comfort and protection



Medent Code	Mft Code	Size	Unit of Sale
13003	73443-04	Padding 7.5cm x 2.4m	Box 12
13002	73443-03	Padding 5.0cm x 2.4m	Box 12
13234	72779-04	Softliner 10.0cm x 10m	Each
13232	72779-02	Softliner 5.0cm x 10m	Each
13233	72779-03	Softliner 7.5cm x 10m	Each



## BSN Dynacast Prelude

### Description

BSN Dynacast Prelude is a synthetic splinting system incorporating polyurethane coated fibre glass covered by polypropylene padding, to offer easy application for the support and immobilisation of fresh fractures and soft tissue injuries. It is presented in a foil covered roll format for flexibility in use and application.

### Features and Benefits

- ✓ A single interlock knit layer of polyurethane resin-impregnated fibreglass
- ✓ The latest innovation in splinting
- ✓ Smoother
- ✓ More conformable
- ✓ Lamination free

Medent Code	Mft Code	Size	Unit of Sale
10314	71415-00	10.0cm x 4.6m	Each
10313	71421-00	7.5cm x 4.6m	Each
10315	71416-00	12.5cm x 4.6m	Each
10317	71419-00	20.0cm x 2.3m	Each
10316	71417-00	15.0cm x 4.6m	Each



## BSN Gypsona\* BP Plaster of Paris Bandages

### Description

Gypsona bandages combine fast setting, high quality plaster of Paris with the conformable and stable interlocked weave of leno cloth. The bandage demonstrates minimum plaster loss while the rich creamy texture allows easy handling.

### Features and Benefits

- ✓ Rapid immobilisation of fractures
- ✓ Leno cloth for good conformability
- ✓ Round cores for improved bandage handling
- ✓ Low plaster loss
- ✓ Moisture proof packaging



Medent Code	Mft Code	Size	Unit of Sale
10402	87060527	Cast 5.0cm x 2.75m	Box 12
10412	87060527	Cast 5.0cm x 2.75m	Each
10406	87061535	Cast 15.0cm x 3.5m	Box 12
10416	87061535	Cast 15.0cm x 3.5m	Each
10403	87067535	Cast 7.5cm x 3.5m	Box 24
10413	87067535	Cast 7.5cm x 3.5m	Each
10404	87061035	Cast 10.0cm x 3.5m	Box 24
10414	87061035	Cast 10.0cm x 3.5m	Each
10346	87004915	Slab 15.0cm x 20m	Each
10344	87004913	Slab 10.0cm x 20m	Each
10347	87004916	Slab 20.0cm x 20m	Each



### 3M Scotchcast Plus Casting Tape

#### Features and Benefits

- ✓ Lightweight for easy mobility, strong and durable casting tape
- ✓ Water-activated resin lets you control working time
- ✓ Allows for weight bearing in 20 minutes
- ✓ Combines the benefits of a fiberglass casting tape with the handling ease of plaster
- ✓ Stretches to conform well to contours
- ✓ Porous weave allows the skin to breathe
- ✓ Contains a lubricating agent in the resin
- ✓ Available in a variety of colours and sizes to maximize patient preferences



Medent Code*	Mft Code	Size	Unit of Sale
10454	82004A	10cm x 3.6m Black	Box 10
10545	82005R	12.5cm x 3.6m Red	Box 10
10555	82005	12.5cm x 3.6m White	Each
10502	82002W	5cm x 3.6m Orange	Box 10
10493	82003G	7.5cm x 3.6m Green	Box 10

\*For illustration purposes only, other colours available for each size.

### 3M Scotchcast Softcast Tape

#### Features and Benefits

- ✓ Fiberglass casting tape that does not set to a completely rigid state, but remains flexible and soft even when it is completely cured
- ✓ Consists of a knitted fiberglass fabric impregnated with a unique polyurethane resin
- ✓ Finished bandage will retain its shape while remaining resilient



Medent Code*	Mft Code	Size	Unit of Sale
10374	82104	10cm x 3.6m White	Pack 10
10372	82102	5cm x 3.6m White	Pack 10
10572	82102B	5cm x 3.6m Blue	Pack 10
10373	82103	7.5cm x 3.6m White	Pack 10

\*For illustration purposes only, other colours available for each size.



## Primenet Polyester Stockinette

### Features and Benefits

- ✓ Cast stockinette made of skin friendly, finest quality synthetic material
- ✓ Use under synthetic cast, Plaster of Paris bandage and protection on amputation stump
- ✓ Protect patients' skin from maceration and irritation during casting
- ✓ Packed in easy to use dispenser carton

Medent Code	Mft Code	Size	Unit of Sale
13312	NET05P	05 Polyester 5cm x 25m	Each
13313	NET07P	7.5 Polyester 7.5cm x 25m	Each
13314	NET10P	10 Polyester 10cm x 25m	Each



## CATEGORY: ORTHOPAEDIC PADDING

### BSN Soffban Orthopaedic Padding

#### Description

Synthetic, protective orthopaedic padding, that blends excellent conformability and good cushioning, with a water-shedding ability that will keep the patients skin dry.

The padding tears quickly and easily around awkward contours, minimising patient trauma and enhancing ease of application for the clinical user

Soffban Natural tears and teases apart easily, aiding rapid application and reducing the risk of stressing the fracture. It is also autoclavable, making it suitable for use in theatre or sterile procedures.

Made entirely from rayon fleece, Soffban# Natural provides soft, comfortable cushioning and protection for the patient, yet still conforms well for easy shaping to all body contours and joints. Ideal for use under Gypsona plaster of Paris bandages, the SoffbanNatural fleece absorbs sweat and exudate, reducing the risk of skin maceration and irritation.

#### Features and Benefits

- ✓ Highly conformable
- ✓ Easily tearable
- ✓ High loft
- ✓ Feathers out
- ✓ Low moisture absorbency



Medent Code	Mft Code	Size	Unit of Sale
13042	71464-01	5CM x 2.7M	Box 12
13044	71486-05	10CM x 2.7M	Box 12
13043	71467-01	7.5CM x 2.7M	Box 12
13046	71465-01	15CM x 2.7M	Box 12
13064	71469-00	Natural 10CM x 2.7M	Box 12
13063	71473-00	Natural 7.5CM x 2.7M	Box 12

### Primepad Polyester Undercast Padding (Soffban Alt)

#### Features and Benefits

- ✓ Undercast padding for cushioning under synthetic cast and Plaster of Paris bandage
- ✓ Made of 100% finest quality polyester
- ✓ Keeps the patient's skin dry and comfortable during long-term wear, easy to tear
- ✓ Smooth surface and good air permeability

Medent Code	Mft Code	Size	Unit of Sale
13022	PAD05P	5cm x 2.7m	Box 12
13023	PAD07P	7.5cm x 2.7m	Box 12
13024	PAD10P	10cm x 2.7m	Box 12
13026	PAD15P	15cm x 2.7m	Box 12



## CATEGORY: CASTING PROTECTOR

### Aquashield Waterproof Cast Protector

#### Features and Benefits

- ✓ Keeps a cast or bandage dry which is critical to the safe healing process
- ✓ Made from heavy-gauge polyurethane, a material renowned for its strength and elasticity
- ✓ Features a pre-formed opening that can be stretched to 5 times its size yet return to within 3% of its original dimension
- ✓ Slides easily over a cast or bandage to create a comfortable & watertight seal



Medent Code	Mft Code	Size	Unit of Sale
25184	A17	Small Half Arm	Each
25185	A21	Regular Half Arm	Each
25204	A26	Small Full Arm	Each
25205	A31	Regular Full Arm	Each
25224	L20	Small Half leg	Each
25225	L25	Regular Half Leg	Each
25244	L33	Small Full Leg	Each
25245	L44	Regular Full Leg	Each

**CATEGORY: BRACE/SPLINT/SUPPORT****Walker Moon Boot****Features and Benefits**

- ✓ Contoured strut design, ensuring a better fit
- ✓ Unique shock absorbing sole reduces impact from the heel strike to improve comfort during ambulation
- ✓ The rocker bottom is specially engineered to be low and wide to help promote a natural gait and increase stability



Medent Code	Mft Code	Size	Unit of Sale
24174	ALL20103830	11 Inch High Small	Each
24175	ALL20103840	11 Inch High Medium	Each
24176	ALL20103850	11 Inch High Large	Each
24204	ALL21101630	17 Inch High Small	Each
24205	ALL21101640	17 Inch High Medium	Each
24206	ALL21101650	17 Inch High Large	Each

**Ankle Stabiliser****Features and Benefits**

- ✓ Helps support & stabilise the ankle
- ✓ Lace up design for a snug and customised fit
- ✓ Non-elastic straps lock the heel to resist inversion & eversion
- ✓ Elastic closure strap 'locks in' tied laces & stabilising straps
- ✓ Fits left or right ankle
- ✓ For mild to moderate ankle sprain, ankle instability, recurrent ankle sprain



Medent Code	Mft Code	Size	Unit of Sale
23574	AOA1830	Small	Each
23575	AOA1830	Medium	Each
23576	AOA1830	Large	Each

**AFO-Deluxe-Allcare****Features and Benefits**

- ✓ Designed to support flacid drop foot
- ✓ Semi rigid calf section for toe clearance and support
- ✓ Thin and flexible foot plate for easy trimming a full length sole



Medent Code	Mft Code	Size	Unit of Sale
23516	ALLAFODLX-LL	L-Large	Each
23546	ALLAFODLX-RL	R-Large	Each
23525	ALLAFODLX-LM	L-Medium	Each
23555	ALLAFODLX-RM	R-Medium	Each
23534	ALLAFODLX-LS	L-Small	Each
23564	ALLAFODLX-RS	R-Small	Each

## BSN Cast Shoe Adult

### Features and Benefits

- ✓ BSN Cast shoe is a premium shoe, constructed of only quality material to aid in mobilisation and to protect the cast from breakdown during periods of weight bearing.



Medent Code	Mft Code	Size	Unit of Sale
25264	58-35010	Small	Each
25265	58-35011	Medium	Each
25266	58-35012	Large	Each

\* Small (5.5 - 7.5 size), Medium (8 - 10 size), Large (9 - 11.5 size)

## Elastic Knee Support/Stabiliser

### Features and Benefits

- ✓ Reduces swelling from minor injuries
- ✓ Flexible compression provides knee support
- ✓ Made of seamless, four-way stretch elastic breathable material for a comfortable fit
- ✓ Open patella
- ✓ For knee pain, articular effusion & swelling, soft tissue injuries, chronic inflammation after operation or trauma



Meden Code	Mft Code	Size	Unit of Sale
23694	AOK1430	Small	Each
23695	AOK1440	Medium	Each
23696	AOK1450	Large	Each

## Arm Sling - Immobiliser with Waist Strap

### Features and Benefits

Has the same features and benefits as that of the sling arm with strap support but has an additional padded cotton waist strap for greater stability.



Medent Code	Mft Code	Size	Unit of Sale
23604	AOS1630	Small	Each
23605	AOS1640	Medium	Each
23606	AOS1650	Large	Each

## Sling Arm with Strap Support

### Features and Benefits

- ✓ Breathable, lightweight and durable material
- ✓ Padded cotton shoulder strap for greater comfort
- ✓ Envelope style with thumb support
- ✓ Universal design fits right or left arm
- ✓ Supports and immobilises the arm following soft tissue injuries, fracture or surgery



Medent Code	Mft Code	Size	Unit of Sale
24024	AOS1730	Small	Each
24025	AOS1740	Medium	Each
24026	AOS1750	Large	Each



## Mesh Arm Sling

### Features and Benefits

- ✓ Supports and immobilises the arm following soft tissue injuries, fracture or surgery
- ✓ Vented nylon mesh with shoulder strap
- ✓ Envelope style with thumb support
- ✓ Universal design fits right or left arm
- ✓ For shoulder, dislocation/subluxation, upper arm fracture, post surgery



Medent Code	Mft Code	Size	Unit of Sale
23914	AOS2330	Small	Each
23915	AOS2340	Medium	Each
23916	AOS2350	Large	Each

## Universal Wrist Splint

### Features and Benefits

- ✓ Metal palmar splint and elastic wrist strap to stabilise the wrist
- ✓ Universal design - fit the left or right hand
- ✓ Provides compression and support to injured and weak wrist
- ✓ For carpal tunnel syndrome, soft tissue injuries, wrist pain, wrist strains, wrist sprains



Medent Code	Mft Code	Size	Unit of Sale
24134	AOW11BE30	Small	Each
24135	AOW11BE40	Medium	Each
24136	AOW11BE50	Large	Each
24137	AOW11BE60	XLarge	Each

## Neoprene Wrist Stabiliser

### Features and Benefits

- ✓ Palmer and dorsal metal splints stabilise the hand and wrist
- ✓ Provides compression and heat retention
- ✓ Supports injured and weak wrist
- ✓ Hook and loop closure for easy adjustment and secure fit
- ✓ Wrist strap with D-ring for additional support and stability
- ✓ Made of 3mm neoprene
- ✓ Available in left or right hand



Medent Code	Mft Code	Size	Unit of Sale
23994	DPNWSS	SL	Each
24004	DPNWSSR	SR	Each
23975	DPNWSM	ML	Each
23985	DPNWSMR	MR	Each
23956	DPNWSL	LL	Each
23966	DPNWSLR	LR	Each

Small 14 - 16cm, Medium Left 16 - 19cm, Large Left 19 - 21.5cm

## Thumb brace

### Features and Benefits

- ✓ Stabilises the joints of thumb, relieves pain
- ✓ Provides thermal support and compression
- ✓ Anatomically contoured metal thumb splint to stabilise the thumb
- ✓ Hook and loop closure for easy adjustment
- ✓ 3mm thick neoprene
- ✓ For joint injury, rheumatoid arthritis, thumb pain, sprain/strain



Medent Code	Mft Code	Size	Unit of Sale
24074	AOW1230	Small	Each
24075	AOW1240	Medium	Each
24076	AOW1250	Large	Each

## Finger Splint Zimmer (Aluminium)

### Features and Benefits

- ✓ For immobilisation of fingers after injury or as a splint during treatment
- ✓ Also used as an arm immobiliser in neonatal treatment



Medent Code	Mft Code	Size	Unit of Sale
23730	11001089	1/2 Inch x 18 Inch	Each
23720	11003089	1 Inch x 18 Inch	Each

## Finger Splints Stax

### Features and Benefits

- ✓ Designed to protect and support the distal joint of the finger
- ✓ Sold individually and come either clear or skin tone
- ✓ Also known as Mallet finger splints



Medent Code	Mft Code	Size	Unit of Sale
23741	11001010006	5cm	Each
23742	11001020006	5.5cm	Each
23743	11001030006	5.8cm	Each
23744	11001040006	6.1cm	Each
23745	11001050006	6.3cm	Each
23748	11001055006	6.8cm	Each
23746	11001060006	7.3cm	Each
23747	11001070006	8.4cm	Each

## Collar Neck Hard Adjustable Adult

### Features and Benefits

- ✓ Ambu Adjustable Neck Collar - Adult (Multifit) - Single



Medent Code	Mft Code	Size	Unit of Sale
23680	COLNECADJB	Adult	Each

## Cervical Collars

### Features and Benefits

- ✓ Soft foam collars designed to limit cervical flexion
- ✓ Made from high density, non latex foam
- ✓ Covered in natural 100% cotton stockinette to provide a comfortable and contoured fit



Medent Code	Mft Code	Size	Unit of Sale
23654	SM400	75 x 480mm. (Small)	Each
23655	SM401	90 x 480mm. (Medium)	Each
23656	SM402	100 x 480mm. (Large)	Each

## Cardboard Splint

### Features and Benefits

The Cardboard Splint is a low cost, compact and convenient splint with a capacity to immobilise a broad range of injuries



Medent Code	Mft Code	Size	Unit of Sale
23636	FRHE-535	Large (L77cm, W33cm)	Each
23638	FRHE-536	Small (L41cm, W25cm)	Each

## CATEGORY: MOBILITY

### Crutches Adult Double Adjustable Forearm (Pair)

#### Features and Benefits

**Ergonomic Crutches.** Ergonomic shaped handle to help distribute the pressure more evenly across the hand. The cuff provides a secure fitting to the forearm and they are adjustable between cuff and handle & handle and ground.

Medent Code	Mft Code	Size	Unit of Sale
25725	CRU691400	Medium (890 X 1200mm)	Pair
25730	CRU691800	Large (890mm -1370mm)	Pair
25705	CRU691500	Ergonomic Med	Pair
25710	CRU691900	Ergonomic Tall	Pair



Ergonomic Crutches

### Stoppers for Crutches (Pair)

- ✓ Standard Crutch Tip
- ✓ Grey non-marking rubber tips
- ✓ High quality rubber tips reinforced with steel to prolong their life
- ✓ Provide good non-slip traction
- ✓ Suit most canes, walking sticks, crutches and anything requiring a non-slipper rubber stopper



Medent Code	Mft Code	Size	Unit of Sale
25750	16MM	Internal 16mm	Pair
25740	CRU686320	Internal 22mm(Single)	Each

## To Order Contact Us at:

### Medent Medical

#### Wellington

Unit 4, 110 Whakatiki St  
PO Box 40356  
Upper Hutt  
New Zealand

**Free Phone: 0800 263 3368**

Phone: 64 4 527 7702

**Free Fax: 0800 527 7499**

Fax: 64 4 5277499

**Email Enquires:** [enquiry@medent.co.nz](mailto:enquiry@medent.co.nz)

#### Auckland

24D Morrin Rd  
Mt Wellington 1072  
Auckland  
New Zealand

**Free Phone: 0800 263 3368**

Phone: 64 9 574 6832

Finding a product and ordering it – you choose an option that suits you.  
Either talk with one of our team or order right here online!

- Order from our sales representative when they visit you.  
To arrange a visit please phone us on 0800 263 3368
- You can phone us on 0800 263 3368  
Or Fax us on 0800 527 7499

All products and more available at [www.medent.co.nz](http://www.medent.co.nz)



---

Notes:

[Online Sales & Ordering](#)



[Email Enquiries](#)

*Medent Barr Guarantee*

If You Are Not 100% Satisfied  
With This BARR'S Branded  
Product Return It for A 100%  
Refund Used or Not

MEDENT  
BARR  
100%

